SENNHEISER

Audioport A 200 Set

Audiology | Stereo Personal Sound Amplifier

Cat. No. 004454

Description

The A 200 is a stereo personal sound amplifier for people who may sometimes need to rely on artificial assistance for their hearing. Due to the stereo amplification via two integrated electret microphones, the user can clearly determine the direction of the sound source. Dynamic headphone transducers ensure a natural and transparent sound reproduction with optimum speech intelligibility. The A 200 features an auxiliary audio input for connecting an external microphone or infra-red receiver.



Features

- Reading glasses for the ears: personal sound amplifier for people with slight to moderate hearing impairment
- Automatic limiter system
- Extremely easy to operate
- Adjustment of frequency response to the listener's hearing
- Practical charger for storing and charging the A 200
- Delivery includes: 1 A 200 personal sound amplifier, 1 BA 151 accupack, 1 L 200 charger for accupack, 1 plug-in mains unit for L 200, 1 spare set of ear-cushions (softer version), 1 set of sealing tape for adjusting the frequency response

Recommended Accessories Spare accupack BA 151	Cat. No. 004146
Replacement Parts Replacement ear-cushions	
(black, firmer version)	Cat. No. 037080
Replacement ear-cushions	.
(transp., softer version)	Cat. No. 040949
Product Variants	
Audioport A 200-120 as Audioport A 200,	
but with L 200-120 charger	Cat. No. 004455
Audioport A 200-UK as Audioport A 200,	
but with L 200-UK charger	Cat. No. 004456

Technical Data

Frequency response	
Transducer principle	dynamic
Microphones	2 electret microphones, stereo
Pick-up pattern	omni-directional
Acoustic amplification	max. 50 dB
Characteristic SPL at 1 kHz	max. 105 dB
Power consumption	5 mA
Operating time	approx. 12 h per accupack
Weight (incl. accupack)	51 g

