

## **Beetle Bluetooth H-1 Headset**

- **Make and receive wireless, hands-free calls via your Bluetooth enabled cell phone**
- **Listen to music and online radio stations via your Bluetooth enabled computer**

**Simply place the T-coil inductive ear hook over your ear, switch your hearing aid or cochlear implant to T-coil, turn on Beetle H-1, and activate Bluetooth on your cell phone or computer.**

- **Wireless, hands-free cell phone communication via Bluetooth and T-coil**
- **Wireless computer CD/DVD audio and VOIP communications such as Skype via Bluetooth and T-coil**
- **Single multi-functional button to answer and end calls**
- **Volume controls customize listening level with impressive amplification**
- **Light weight, inductive ear hooks provide strong signal strength with reduced background noise**
- **T-coil/Telecoil compatible - BTE, ITE and ITC aids**

### **Beetle H-1 specifications**

- **Bluetooth specification: v1.2**
- **Bluetooth profiles: Headset (HSP) and handsfree (HFP)**
- **Weight: 18 grams**
- **Dimension: 50mm(L), 24mm(W),19mm(H)**
- **Talk time: Up to 6 hours \***
- **Stand-by time: Up to 280 hours \***
- **Operating range: Up to 10 meters**
- **Radio Frequency: 2.4-2.4835 GHz ISM Band**
- **Transmitter Power: Class 2, ~1mW**
- **Power Consumption: 27mA (working); <1mA (stand-by)**
- **Battery: Rechargeable Lithium-ion polymer battery 3.7V**
- **Charger: USB cable**
- **Charging Time: 1.5 hours - 2 hours**
- **Certification: FCC, CE, RoHS, BQB**
- **1-year manufacturer warranty**

**\* May vary due to use, network and mobile phone settings**