



conversOR

The perfect companion for your hearing instrument

- Greatly reduces unwanted background noise and enhances sound level and clarity.
- Enjoy easy, relaxed conversations – The sound you hear is clearer and annoying interference is minimized.
- With Zoom and Wide settings you can focus on single or multiple sound sources (for one-to-one or group conversations)
- Use it anywhere – In the car; in classroom settings, connect it to your TV, take it to the theater, or to any social gathering.
- Easy, discreet use with a variety of Telecoil-equipped hearing instruments: BTE (Behind the Ear), ITE (In the Ear), ITC (In the Canal) and some Cochlear implants.
- Microphone to Receiver transmission up to 150 feet!



Lightweight Pendant Receiver


Wireless FM Microphone/Transmitter

listenOR

Enhanced sound quality for those without a hearing instrument*

- Amplifies available sound level and greatly reduces unwanted background noise.
- Provides you greater hearing pleasure, even in noisy environments.
- Ideal for use with TV or other audio devices, or when conversing with friends and family.
- Zoom and Wide settings allow you to focus either on single or multiple sound sources.
- Includes a pair of earphones with volume control.

* Can also be used with Telecoil-equipped hearing instruments with an optional loop, ordered separately.



Powerful Microphone with Earphones