

A smiling man and woman are shown in a close-up shot. The man is in the foreground, wearing a white sweater and a green scarf, looking towards the camera with a wide smile. The woman is behind him, wearing a white sweater and a black fur vest, also smiling. The background is a soft-focus field of golden autumn leaves. In the top right corner, there is a logo for Williams Sound with the tagline 'helping people hear'. The bottom right of the image features a light blue geometric graphic with the text 'assistive listening solutions' and the slogan 'improve your hearing, improve your life.'

WILLIAMS SOUND®
helping people hear

assistive
listening solutions

improve your hearing, improve your life.



[“helping people hear” is not just our motto, it is the core of our business philosophy. We understand the impact and importance of hearing in people’s everyday lives - that’s why we set our sights to be a legendary company with remarkable products and exceptional service.]

assistive listening devices

Assistive listening devices make it easier to hear in everyday situations, helping you “focus” on the sounds and voices most important to you. They reduce background noise and compensate for poor room acoustics or distance from the sound source — wherever you go, whatever you do. In restaurants, classrooms or personal conversations, you’ll hear easily and clearly as sound and voice travel directly to your ear. Listening is significantly improved on the phone or while watching TV. Assistive listening devices work with mild to severe hearing loss — with or without a hearing aid. So, relax. Enjoy the sounds of life and stay connected to your world.

To learn more about assistive listening solutions visit us at www.williamssound.com. Call us toll free at 800.843.3544.



improve your hearing improve your life.

Improve your hearing with the Pocketalker® personal amplifier. The Pocketalker® amplifies sounds closest to the listener while reducing background noise. Ideal for one-on-one conversation, small group and television listening or conversing in the car. Use with or without hearing aids. Simply plug in your earpiece, position the mic near the preferred sound, adjust volume to your comfort — and start listening!

- 200 hours of battery life
- Adjustable tone and volume control (on ULTRA model)
- Choice of earphones or headphones
- Use with neckloops and telecoil-equipped hearing aids
- Amplifies sound for better understanding
- 5 year warranty — 90 days on accessories



Pocketalker® Ultra
model PKT D1

personal amplifiers



Pocketalker® Pro
model PKT PRO



reduce background noise maximize listening.

Motiva™ system delivers superior audio directly to the ear, enhancing sound at home, work or school - any place where background noise or distance from the source makes listening difficult. Ideal for classrooms, group listening and training. Simply plug your earpiece into the receiver, position the transmitter near the sound source, adjust volume to your comfort — and start listening!

- PFM 360 receiver features built-in Pocketalker®
- Portable, RF wireless operation
- Overcomes distance from speaker to listener
- 150 foot range
- Operates up to 60 hours at a time
- 5 year warranty — 90 days accessories
- Rechargeable systems available

personal fm listening



Motiva® Personal Listening System
model PFM 360



Motiva® Personal Listening System
model PFM 330



your clear connection to bluetooth technology.

CM-BT provides everything you need to listen to your mobile phone, computer or other Bluetooth-enabled device with your telecoil-equipped hearing aids. With a push of a button, calls via a mobile phone on a belt clip, in a pocket or in a purse can be answered / terminated.

- Compatible with all three Bluetooth protocols
- Provides a 40 dB low-to-high volume range when used with a neckloop
- Provides up to 118 dB SPL when used with a headphone or earphone
- Range of up to 35 feet
- Provides up to 6 hours of talk time
- Supports Bluetooth specification 1.2
- Covered by a 2-year warranty

bluetooth amplifier



CM-BT Neckloop Amplifier
model CM-BT



the “marriage saver” listen without disturbing others.

Stylish and lightweight, the TVTalker™ listening system is easy to use and affordable. Use the headset to listen to the TV or stereo without disturbing others. Discover why they call this the “marriage saver.”

- Stereo operation, 2.3/2.8 MHz
- Lightweight, under-the-chin style receiver (WIR 240)
- Plugs directly into your TV, DVD player, or stereo
- Receiver charges itself when placed back on top of transmitter
- Automatic shut off (WIR 240)
- 2 year warranty — 90 days on accessories
- FM model (WFM 260) now available. Freedom to listen in 100' range



TV Talker™
WIR 240

TV listening



TV Talker™ Flex
WIR 250

www.tvtalker.net



listen clearly never miss a call.

The Teletalker™ amplified telephone features the strongest amplification available on the market, 55 dB. It improves speech understanding and amplifies calls without feedback.

- 55 dB of strong amplification
- Improves speech understanding
- Extra loud ringer, flashing ring indicator
- Easy-to-use buttons and volume control
- Compatible with telecoil-equipped hearing aids
- 5 year warranty — 90 days on accessories



amplified telephones

The Ringmax is an amplified phone ringer that connects to your phone and lets you customize your phone ring for better hearing. Quickly and easily increase the volume, then set the tone and pitch that's best for you.

- 95 dB of amplification
- Four selectable tones / warble
- Includes flashing ring indicator
- 2 year warranty — 90 days on accessories





EAR 008



EAR 013



EAR 014



EAR 041



EAR 022



HED 021



HED 026



HED 027



MIC 014-R



MIC 027



MIC 049



MIC 054



MIC 090



MIC 100



NKL 001



BAT KT3



BAT KT6



BAT AP07D

accessories

EAR 008

Wide Range Earphone. Hangs on outside of ear. Excellent sound quality. 32 Ω , mono

EAR 013

Single Mini Earbud. Fits in the outer part of the ear. 32 Ω , mono

EAR 014

Dual Mini Earbud. Fits in the outer part of the ear. 16 Ω , mono

EAR 041

New Single Mini Isolation Earbud. Fits in the outer part of the ear. 32 Ω , mono

EAR 022

Surround Earphone. Hangs on outside of ear. Excellent sound quality. 32 Ω , mono

HED 021

Folding headphone. Adult size. 32 Ω , mono

HED 026

Deluxe rear-wear headphone. Adult size. 16 Ω , mono

HED 027

Deluxe folding headphone. Adult size. 16 Ω , mono

MIC 014-R

Plug mount mic. No cord. Omnidirectional Condenser

MIC 027

Hand-held unidirectional mic. Improved RF resistance. Unidirectional Condenser, 20-18 kHz

MIC 049

Omnidirectional conference mic. Enhance pickup range (8-12ft), 36" cord. Omnidirectional Condenser, 100-8 kHz

MIC 054

Directional lapel clip mic, mono. Unidirectional Condenser

MIC 090

Mini lapel clip mic. 3.5mm plug, 39" cord, mono. Omnidirectional Condenser, 20-18 kHz

MIC 100

Unidirectional headband mic. Improved RF resistance. Ideal for high noise environments. Unidirectional Condenser

NKL 001

Neckloop telecoil coupler. For use with hearing aids equipped with a T-coil switch or an induction earphone. 3.5 mm plug, mono, 8 - 16 Ω

BAT KT3

3V charger unit and (2) AA NiMH batteries (BAT 026)

BAT KT6

3V charger unit (CHG 3502) and (4) AA NiMH batteries (BAT 026)

BAT APO7D

Rechargeable battery for TV Talker™ receivers, NiMH



For more information about Williams Sound's full line of
assistive listening products and accessories, contact our sales department at
800.843.3544 / info@williamssound.com / www.williamssound.com