

Comfort Contego®

A digital wireless assistive listening device



Bringing you closer together

Comfort Contego® – a digital wireless assistive listening device improves the quality of life of people with hearing loss in many everyday situations.

The clear sound and unique features of Comfort Contego® have contributed to its rapid success in Europe and the US. Comfort Contego® sends its signals encrypted to prevent electronic eavesdropping, an industry first by Comfort Audio.

Both the transmitter and the receiver have built-in zoom microphones that can be adjusted by the user for maximum sound customization. Depending on conditions, the microphones can be set for omni or directional use, and volume levels can be adjusted remotely.



Features

-  In meetings at work
-  In the classroom
-  Watching TV
-  Talking on the phone
-  Secure personal privacy



Comfort Audio
HEAR THE FUTURE

In everyday situations




At a lecture



Talking on the phone



At a meeting



Watching TV

Technical facts

Transmission: Digital modulation with encryption.

Range: Seventy five (75) feet.

Size: 0.87 x 1.77 x 3.86 inches

Weight: Weight: 2.1 oz.

Antenna: Built in.

Channels: 38

Battery: Rechargeable Li-ion.

Charging time: Four (4) hours.

Operating time: Sixteen (16) hours.

Microphones: Built in - directional and omni-directional.

Frequency range: 904-926 MHz.

Input: Line in/Microphone.

Max Output: <130 dB SPL

Max Gain: <50 dB



Comfort Audio
HEAR THE FUTURE