

Westone® Custom-Fit Earpiece Order Form

1.800.525.5071 / westone.com

Please Note: Fill out order form completely; failure to do so may delay the order

Patient Information

Name _____ Date ____/____/____

Dispenser Information

Bill To: The invoice for this order will be sent to this address

Submitted By _____

Westone Acct No. _____

Shipping Code _____

P.O. No. _____

Medicaid Information

Medicaid No. _____

Billing Address

Address _____ Phone _____

City _____ State _____ Zip _____

Email _____

Shipping Instructions

Ship To: (use address listed above)

Dispenser Other Office Patient

Name _____

Address _____

City _____ State _____ Zip _____

Shipping Method

FedEx® Priority Overnight® (10:30 a.m. next day) FedEx Standard Overnight® (3:00 p.m. next day)

FedEx 2Day® (4:30 p.m. second business day) FedEx Express Saver® (4:30 p.m. third business day)

Standard Shipping (1st Class Mail®)

Rush Service

Yes, rush this order (\$15.00 rush charge will apply to each earpiece only, outbound charges are not included)

Billing Method

Westone Account Check or Money Order Enclosed Credit Card

Credit Card Information

Billing Address _____

City _____ State _____ Zip _____

Phone _____ Ext. _____

Card No. _____ Exp. ____/____/____

CVV Code _____

Printed Name _____

Signature _____

(Required if using credit card as payment method)

Special Requests / Information

(If a remake – please provide reason)



Right Left (Use back of order form for additional notes)

Yes, Westone may make necessary changes to an earpiece style or material without contacting me.

Please Send More: Order Forms Impression Boxes

Track this order online: webear.westone.com

Earpiece Style/Information

Style No. Left _____ Right _____

Earpiece

Helix Lock Linda Helix
 3/4 Shell Standard Earhook
 RIC Fitting Slim-Tube

Manufacturer _____

Model _____ Power (dB) _____

Canal Length

Short Medium Long
 Do Not Cut As Marked on Imp

Earpiece Materials and Colors

Visit: westone.com for complete color chart

SILICONE

OtoBlast™ W-1™
 Single Color Swyrl™ (3 colors)
Color 1 _____ Color 2 _____
Color 3 _____

DisappEar™ "A" "B" "C"
(Please Note: "C" not available in W-1™)

Cat-Eyes

Clear + 2 Colors
Color 1 _____ Color 2 _____

Glitter

Single Color Multi-Color
 Mix (3 colors) Color 1 _____
Color 2 _____ Color 3 _____

Canal Additive

e-Compound™ (Matte finish only)

Earpiece Finish

Gloss Matte Slick Sil®

AQUANOT™ (Not available in Clear)

Single Color Swyrl™ (3 colors)
Color 1 _____ Color 2 _____
Color 3 _____

Earpiece Finish

Gloss Matte Slick Sil®

ACRYLIC

Single Color Swyrl™ (3 colors)
Color 1 _____ Color 2 _____
Color 3 _____

DisappEar™ "A" "B" "C"

Canal Additive

Flex Canal (Beige/Soft)
 e-Compound™ (White/Hard)
 Flex-com (Tan/Semi-Soft)

Earpiece Finish

Gloss Matte

Glitter

Single Color Multi-Color
 Mix (3 colors) Color 1 _____
Color 2 _____ Color 3 _____

VINYL (PVC)

Formula II

Light Tan Light Brown
 Medium Brown Dark Brown

Formula II Clear

Clear Color _____


Earpiece Finish

Gloss Matte

Patient Hearing Loss

Mild Moderate
 Severe Profound

(Please Complete Audiogram)



Left Ear – X Right Ear – O

Venting Options (Select One)

Vent Size	Metric	Decimal
<input type="checkbox"/> Pressure	1.0 mm	0.040"
<input type="checkbox"/> Standard	1.6 mm	0.063"
<input type="checkbox"/> Medium	2.0 mm	0.081"
<input type="checkbox"/> Large	2.4 mm	0.096"

MVP (includes all sizes)
 SAV* (includes all sizes)
 Largest vent possible
 Semi-IROS
 No Vent

Please Note: If a vent is not selected, earpiece will be made with no vent. RIC earpieces will receive an MVP if no vent is selected.

Semi-IROS – Vent side of canal tip is cut back by approximately 1/2.

*SAV – Used only if space allows. MVP may be used instead.

Tubing Options

Factory Select **or** Tube # _____
 DisappEar™ Tubing
 Dry Tube
 Tube Through (if size permits)
 Permanently Glue Tube
 TRS **or** Tube Lock
 Friction Fit Tube
 No Tube

Additional Information

Initials (3 max.)

Patient Initials _____

Right and Left Side Indicators

Red (R) / Blue (L) Initials
 Red (R) / Blue (L) Color Dots

Removal Handles

Post Type Handle(s)
 Removal Filament(s)
 Half Moon Handle

Earpiece Cords

Nylon Glue Cord – Black
 Nylon Pop Cord – Black
(Please Note: Black is the default color)

Metal Detection

Embedded Metal Required

DefendEar® Packaging

Please Use: DefendEar® Industrial Packaging for this order.

Fold Here

Fold Here